



End-to-end solution for a

quick functional assessment

In less than 3 minutes Lynx evaluates the functional capacity of the patient

IN JUST
3 STEPS



MEASURE

Place the Bluetooth sensors on the joint to monitor the patient's movements

2

ANALYZE

With a simple click, Lynx automatically quantifies all the key parameters of the movement.

3



INFORM

Lynx generates a complete patient report with graphical and numerical representation of the results.

What are the advantages of Lynx?



- 3D analysis of motion limitations
- Accurate and real-time biofeedback
- Simultaneous measurement of 2 bilateral joints
- Simple and intuitive results
- Quality assurance: CE marking of medical device



Lynx is simple, fast and accurate

COMPARISON MODULE

An objective comparison of the functional capacity of the patient helps in clinical decision making.

EVOLUTION MODULE

The tests performed on patients are stored to monitor their evolution over time, thus facilitating the follow-up of their progress.

Discover Lynx products



Lynx DYNAMO

Strength analysis module.
Using a single sensor, tension forces of up to 100 kg can be measured.

- Maximum and average load
- Response time at maximum load
- Holding time at maximum load

Lynx ROM

Joint motion analysis module. The sensors can be attached to any joint in the body.

- Maximum and average Range of Motion (ROM)
- Speed
- Movement Harmony
- Position vs Velocity graph



















